

Anti-human TAU total mab 12C2

Catalogue number	Package size
0102005201	100 µg
0102005203	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	12C2
Host:	mouse
Isotype:	IgG2a _k
Immunogen:	recombinant human TAU 441
Specificity:	human TAU
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at –20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	PHF TAU, TAU 441

