

Anti-human phospho-181 TAU mab 1E7

Catalogue number	Package size
847-0102003801	100 µg
847-0102003803	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	1E7
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Synthetic peptide CIPAKTPPAPKT(Pi)PPS
Specificity:	human TAU phosphorylated at amino acid position 181 (T)
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	PHF TAU, phosphorylated TAU at amino acid position 181 (T)

