

Anti-human TAU-E3 mab 2B6

Catalogue number	Package size
847-0102006401	100 µg
847-0102006403	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	2B6
Host:	mouse
Isotype:	IgG1
Immunogen:	Recombinant human TAU 441 in vitro aggregated
Specificity:	Exon 3 of human TAU441, amino acids 74-102
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	Phosphate buffered saline ph 7.4, without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot
Dilution:	1 µg/ml in Western Blot
Reactivity:	Tau441, Tau410

