

Anti-human prion protein mab 6E2

Catalogue number	Package size
847-0102004101	100 µg
847-0102004103	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	6E2
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	recombinant human prion protein
Specificity:	human prion protein
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot (see figure below)
Dilution:	1-10 ng/ml in ELISA, 2 µg/ml in Western Blot
Reactivity:	prion protein from human and cattle

