

Anti-human β -amyloid(42) mab 10E2**Catalogue number**847-0102012001
847-0102012003**Package size**100 μ g
1 mg*For research use only.*

Product:	monoclonal antibody
Clone:	10E2
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Synthetic peptide SGGVVIA
Specificity:	β -amyloid-42, specific to C-terminal with amino acids 41+42
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	2 ng/ml in ELISA, 1 μ g/ml in Western Blot
Reactivity:	to C-terminal with amino acids 41+42

Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • GermanyCEO
Dr. Ingolf LachmannPhone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.comUst-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig