
Anti-human phosphoTDP43 Ser409/Ser410 **mab 11C10**

Catalogue number	Package size
847-0102008301	100 µg
847-0102008303	1 mg

For research use only

Product:	monoclonal antibody
Clone:	11C10
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	TDP43-Peptide phosphorylated at Ser409 and Ser410
Specificity:	phosphoTDP43 Ser409/Ser410
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, IHC
Dilution:	0.1 - 1 µg/ml in Western Blot and IHC
Reactivity:	phosphoTDP43



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

CEO
Dr. Ingolf Lachmann

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig