
Anti-human TAU MTBR mab 12C8**Catalogue number**

847-0102009801

847-0102009803

Package size

100 µg

1 mg

For research use only.

Product:	monoclonal antibody
Clone:	12C8
Host:	mouse
Isotype:	IgG2b _k
Immunogen:	human TAU sequence VYKPVDLSKVTSCG
Specificity:	VYKPVDLSKVTSCG human TAU (2N4R) amino acids 309-323
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	Carbonate buffer pH 9.6 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	Microtubuli binding region TAU protein



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

CEO
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig