

Anti-human phospho-181 TAU mab 14B6

| Catalogue number | Package size |
|-------------------------|---------------------|
| 847-0102008601 | 100 µg |
| 847-0102008603 | 1 mg |

For research use only.

| | |
|----------------------|---|
| Product: | monoclonal antibody |
| Clone: | 14B6 |
| Host: | mouse |
| Isotype: | IgG1 _k |
| Immunogen: | Synthetic peptide CIPAKTPAPKT(Pi)PPS |
| Specificity: | human TAU phosphorylated at amino acid position 181 (T) |
| Purification: | purified by affinity chromatography |
| Purity: | > 95% by SDS-PAGE |
| Solution: | PBS pH 7.4 without additives |
| Storage: | at -20 °C (repeated thawing and freezing should be avoided) |
| Applications: | ELISA, Western Blot, Immunohistochemistry |
| Dilution: | 1 µg/ml in Western Blot |
| Reactivity: | PHF TAU, phosphorylated TAU at amino acid position 181 (T) |

