

Anti-human phospho-199 TAU mab 1C7

Catalogue number	Package size
0102011601	100 µg
0102011603	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	1C7
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Synthetic peptide. <i>Design and synthesis were kindly performed by group of Prof. Dr. Ralf Hoffmann (BBZ, University of Leipzig, Germany).</i>
Specificity:	human TAU sequence SGYSSPGSPGTPG phosphorylated at amino acid position 199
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	PHF TAU, phosphorylated TAU at amino acid position 199



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

CEO
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig