

Anti-human phospho-235 TAU mab 1F5

Catalogue number	Package size
0102011901	100 µg
0102011903	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	1F5
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Synthetic peptide KKVAVVRT(Pi)PPKSPSS
Specificity:	human TAU phosphorylated at amino acid position 235 (S)
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at –20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	PHF TAU, phosphorylated TAU at amino acid position 235 (S)

