

***Anti-human TAU MTBR mab 1F6***

| Catalogue number | Package size |
|------------------|--------------|
| 847-0102010701   | 100 µg       |
| 847-0102010703   | 1 mg         |

*For research use only.*

---

|                      |                                                             |
|----------------------|-------------------------------------------------------------|
| <b>Product:</b>      | monoclonal antibody                                         |
| <b>Clone:</b>        | 1F6                                                         |
| <b>Host:</b>         | mouse                                                       |
| <b>Isotype:</b>      | IgG1 <sub>k</sub>                                           |
| <b>Immunogen:</b>    | human TAU sequence K18CF4R                                  |
| <b>Specificity:</b>  | human TAU (2N4R) amino acids 244-372                        |
| <b>Purification:</b> | purified by affinity chromatography                         |
| <b>Purity:</b>       | > 95% by SDS-PAGE                                           |
| <b>Solution:</b>     | Carbonate buffer pH 9.6 without additives                   |
| <b>Storage:</b>      | at -20 °C (repeated thawing and freezing should be avoided) |
| <b>Applications:</b> | ELISA, Western Blot, Immunohistochemistry                   |
| <b>Dilution:</b>     | 1 µg/ml in Western Blot                                     |
| <b>Reactivity:</b>   | Microtubuli binding region TAU protein                      |

---



Roboscreen GmbH  
Hohmannstraße 7  
04129 Leipzig • Germany

Phone: +49 341 989 734 0  
Fax: +49 341 989 734 199  
[info@roboscreen.com](mailto:info@roboscreen.com)  
[www.roboscreen.com](http://www.roboscreen.com)

CEO  
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314  
County Court Leipzig  
HRB 17503  
Leipzig