

***Anti-human TAU total mab 1G2***

<b>Catalogue number</b>	<b>Package size</b>
847-0102007501	100 µg
847-0102007503	1 mg

*For research use only.*

---

<b>Product:</b>	monoclonal antibody
<b>Clone:</b>	1G2
<b>Host:</b>	mouse
<b>Isotype:</b>	<i>IgG2a<sub>k</sub></i>
<b>Immunogen:</b>	<i>Aggregated and phosphorylated human TAU 441</i>
<b>Specificity:</b>	<i>TPP motif at position T175 and T181 of TAU441 w/o phosphor</i>
<b>Purification:</b>	purified by affinity chromatography
<b>Purity:</b>	> 95% by SDS-PAGE
<b>Solution:</b>	PBS buffer pH 7.4 without additives
<b>Storage:</b>	at -20 °C (repeated thawing and freezing should be avoided)
<b>Applications:</b>	ELISA, Western Blot, Immunohistochemistry
<b>Dilution:</b>	1 µg/ml in Western Blot
<b>Reactivity:</b>	<i>Tau w/o phosphor at T175 and T181, all isoforms of human TAU</i>

---

