

Anti-human TAU total mab 2B8

Catalogue number	Package size
847-0102007601	100 µg
847-0102007603	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	2B8
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	human TAU 441
Specificity:	Epitope: EAAGHTVQ, human TAU amino acids 117-124
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS ph 7.4 w/o additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot
Reactivity:	TAU 441 all isoforms, TAU total in CSF

