


## Anti-human Brain derived TAU mab 2B8

Catalogue number	Package size
847-0102007601	100 µg
847-0102007603	1 mg

*For research use only.*

<b>Product:</b>	monoclonal antibody
<b>Clone:</b>	2B8
<b>Host:</b>	mouse
<b>Isotype:</b>	IgG1 <sub>k</sub>
<b>Immunogen:</b>	human TAU441 (2N4R)
<b>Specificity:</b>	Epitope: HVTQARMV, human TAU441 (2N4R) amino acids 121-128



<b>Purification:</b>	purified by affinity chromatography
<b>Purity:</b>	> 95% by SDS-PAGE
<b>Solution:</b>	PBS ph 7.4 w/o additives
<b>Storage:</b>	at -20 °C (repeated thawing and freezing should be avoided)
<b>Applications:</b>	ELISA, Western Blot, Immunohistochemistry, Immunoprecipitation, SIMOA, Luminex
<b>Reactivity:</b>	TAU441 (2N4R) all isoforms, TAU total in CSF, brain derived TAU

