

Anti-human TDP43 mab 2C9

Catalogue number	Package size
847-0102010101	100 µg
847-0102010103	1 mg

For research use only

Product:	monoclonal antibody
Clone:	2C9
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	recombinant human TDP43
Specificity:	Amino acids KDNKRKMDETDA of TDP43 protein
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, IHC
Dilution:	0.1 - 1 µg/ml in Western Blot and IHC
Reactivity:	TDP43



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

CEO
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig