

Anti-human β -amyloid mab 3C4

Catalogue number	Package size
847-0102008501	100 μ g
847-0102008503	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	3C4
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Synthetic peptide VKMDAEFRHDSGY
Specificity:	β -amyloid 1-38, 1-40, 1-42
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	2 ng/ml in ELISA, 1 μ g/ml in Western Blot
Reactivity:	β -amyloid



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

CEO
Dr. Ingolf Lachmann

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig