

Immuno Beads Human Beta Amyloid

Catalogue number

847-0601000401
847-0601000402

Package size

6 x IP
24 x IP

For research use only.

Product:	Magnetic beads coupled antibody
Isotype:	Mouse IgG
Specificity:	N-terminus of human beta amyloid peptides
Solution:	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
Concentration:	10 mg magnetic beads/ml and 200 µg IgG/ml
Storage:	2-10°C <i>After delivery of the beads, make sure that the beads are completely in suspension and do not stuck to the tube walls or in the lid. To do this, vortex the beads and shortly spin them down by centrifuge so that all beads and storage liquid are at the bottom of the tube.</i>
Use for:	Immunoprecipitation of human beta amyloid peptides
Application:	25 µl Immuno Beads per sample
Use with:	NEURO-IP Kit (847-0108000108)

