

*Immuno Beads Human  $\alpha$ -Synuclein total***Catalogue number**

847-0601000601

847-0601000602

**Package size**

6 x IP

24 x IP

*For research use only.*

---

<b>Product:</b>	Magnetic beads coupled antibody
<b>Isotype:</b>	Mouse IgG
<b>Specificity:</b>	Human $\alpha$ -Synuclein total 140 amino acids
<b>Solution:</b>	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
<b>Concentration:</b>	10 mg magnetic beads/ml and 200 $\mu$ g IgG/ml
<b>Storage:</b>	2-10°C <i>After delivery of the beads, make sure that the beads are completely in suspension and do not stuck to the tube walls or in the lid. To do this, vortex the beads and shortly spin them down by centrifuge so that all beads and storage liquid are at the bottom of the tube.</i>
<b>Use for:</b>	Immunoprecipitation of human $\alpha$ -Synuclein oligomers
<b>Application:</b>	25 $\mu$ l Immuno Beads per sample
<b>Use with:</b>	NEURO-IP Kit (847-0108000108)

---



Roboscreen GmbH  
Hohmannstraße 7  
04129 Leipzig • Germany

Phone: +49 341 989 734 0  
Fax: +49 341 989 734 199  
[info@roboscreen.com](mailto:info@roboscreen.com)  
[www.roboscreen.com](http://www.roboscreen.com)

CEO  
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314  
County Court Leipzig  
HRB 17503  
Leipzig