

Immuno Beads Human brain derived Tau

Catalogue number	Package size
847-0601000301	6 x IP
847-0601000302	24 x IP

For research use only.

Product:	Magnetic beads coupled antibody
Isotype:	Mouse IgG
Specificity:	Human brain derived Tau, TAU441 (2N4R)
Solution:	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
Concentration:	10 mg magnetic beads/ml and 200 µg IgG/ml
Storage:	2-10°C <i>After delivery of the beads, make sure that the beads are completely in suspension and do not stuck to the tube walls or in the lid. To do this, vortex the beads and shortly spin them down by centrifuge so that all beads and storage liquid are at the bottom of the tube.</i>
Use for:	Immunoprecipitation of human brain derived Tau protein
Application:	25 µl Immuno Beads per sample
Use with:	NEURO-IP Kit (847-0108000108)

