

Immuno Beads Human patho-oligomeric α -Synuclein

Catalogue number

847-0601000501

847-0601000502

Package size

6 x IP

24 x IP

For research use only.

Product:	Magnetic beads coupled antibody
Isotype:	Mouse IgG
Specificity:	Human α -Synuclein oligomers only, amino acids 47-52 (EGVVGVA) in β -sheet
Solution:	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
Concentration:	10 mg magnetic beads/ml and 200 μ g IgG/ml
Storage:	2-10°C <i>After delivery of the beads, make sure that the beads are completely in suspension and do not stuck to the tube walls or in the lid. To do this, vortex the beads and shortly spin them down by centrifuge so that all beads and storage liquid are at the bottom of the tube.</i>
Use for:	Immunoprecipitation of human α -Synuclein oligomers
Application:	25 μ l Immuno Beads per sample
Use with:	NEURO-IP Kit (847-0108000108)

