

Immuno Beads Human total Tau

| Catalogue number | Package size |
|-------------------------|---------------------|
| 847-0601000101 | 6 x IP |
| 847-0601000102 | 24 x IP |

For research use only.

| | |
|-----------------------|---|
| Product: | Magnetic beads coupled antibody |
| Isotype: | Mouse IgG |
| Specificity: | total tau protein all isoforms, big tau |
| Solution: | PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %) |
| Concentration: | 10 mg magnetic beads/ml and 200 µg IgG/ml |
| Storage: | 2-10°C <i>After delivery of the beads, make sure that the beads are completely in suspension and do not stuck to the tube walls or in the lid. To do this, vortex the beads and shortly spin them down by centrifuge so that all beads and storage liquid are at the bottom of the tube.</i> |
| Use for: | Immunoprecipitation of human total tau protein |
| Application: | 25 µl Immuno Beads per sample |
| Use with: | NEURO-IP Kit (847-0108000108) |

