

### *Immuno Beads brain derived Tau*

<b>Catalogue number</b>	<b>Package size</b>
847-0601000301	6 x IP
847-0601000302	24 x IP

*For research use only.*

---

<b>Product:</b>	Magnetic beads coupled antibody
<b>Isotype:</b>	Mouse IgG
<b>Specificity:</b>	Brain derived Tau, TAU441 (2N4R)
<b>Solution:</b>	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
<b>Concentration:</b>	10 mg magnetic beads/ml and 200 µg IgG/ml
<b>Storage:</b>	2-10°C
<b>Use for:</b>	Immunoprecipitation of brain derived Tau protein
<b>Application:</b>	25 µl Immuno Beads per sample
<b>Use with:</b>	NEURO-IP Kit (847-0108000108)

---

