

Immuno Beads p50-Tau

Catalogue number	Package size
847-0601000201	6 x IP
847-0601000202	24 x IP

For research use only.

Product:	Magnetic beads coupled antibody
Isotype:	Mouse IgG
Specificity:	Tau protein phosphorylated at amino acid T50
Solution:	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
Concentration:	10 mg magnetic beads/ml and 200 µg IgG/ml
Storage:	2-10°C
Use for:	Immunoprecipitation of p50-Tau protein
Application:	25 µl Immuno Beads per sample
Use with:	NEURO-IP Kit (847-0108000108)



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

CEO
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig