

Product Launch of the Rapid INSTANT Virus RNA/DNA Kit - FX

Hereby, please be informed, that Roboscreen GmbH launches its new product Rapid INSTANT Virus RNA/DNA Kit - FX for automated nucleic acid extraction with CyBio FeliX by November 1st, 2022.

Rapid INSTANT Virus RNA/DNA Kit – FX



- CE-IvD labeled extraction of viral RNA and DNA
- Automated MAG bead based
- Up to 96 samples per run in appr. 35 min

Intended Use

The Rapid INSTANT Virus RNA/DNA Kit - FX is intended for automated purification of viral nucleic acids from cell free fluid biological samples. The Rapid INSTANT Virus RNA/DNA Kit - FX is designed for extraction using the CyBio FeliX Basic Unit in combination with the CyBio FeliX Extraction Set. With a starting volume of 400 µl per specimen a throughput of up to 96 samples is possible at once. The Rapid INSTANT Virus RNA/DNA Kit - FX is optimized for minimal processing time.

Product Specifications

Basis	Magnetic particle based
Handling	Automated extraction using the CyBio FeliX
Starting material	<ul style="list-style-type: none"> ▪ 400 µl sample volume ▪ Fresh or frozen cell-free biological fluids (e.g. serum, plasma, sputum, swab supernatant)
Duration of extraction	~ 35 min (96 samples)

Order Information

Order Number	Product description	Kit size
847-0259200906	Rapid INSTANT Virus RNA/DNA Kit - FX	2x 96 reactions for 400 µl sample volume
Necessary laboratory equipment		
847-0259200924	Prefilling Set 400 - FX	96 reactions
847-0259200936	Plate Set 400 Rapid - FX	96 reactions
Disposables		
847-FX-TIPS-1000	FX Filter Tips 1000 µl	16x 96 tips

Associated Documents

- [Instruction for use](#)
- [Catalogue Molecular Diagnostic Products](#)
- [Declaration of Conformity](#)

For further information on IvD certified kits for nucleic acid extraction, detection of viral pathogens and products for markers of neurodegenerative diseases please visit our homepage www.roboscreen.com.