

Anti-Fluoresceinisothiocyanat (FITC) mab 7F4

Catalogue number

847-0102006001

847-0102006003

Package size

100 µg

1 mg

For research use only.

Product:	monoclonal antibody
Clone:	7F4
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Conjugate of bovine serum albumin and Fluoresceinisothiocyanat (FITC)
Specificity:	FITC
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot
Dilution:	1 µg/ml in Western Blot
Reactivity:	FITC, PCR products labelled with FITC, proteins labelled with FITC

