

Anti-HANTA mab 5D8

Catalogue number	Package size
847-0102005901	100 µg
847-0102005903	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	5D8
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Recombinant HANTA virus protein
Specificity:	HANTA Dobrava and Puumala strains
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot
Dilution:	1 µg/ml in Western Blot
Reactivity:	HANTA virus protein strains Dobrava and Puumala

