

Anti-human β -amyloid(42) mab 22C2

| Catalogue number | Package size |
|------------------|--------------|
| 847-0102008901 | 100 μ g |
| 847-0102008903 | 1 mg |

For research use only.

| | |
|----------------------|--|
| Product: | monoclonal antibody |
| Clone: | 22C2 |
| Host: | mouse |
| Isotype: | IgG1 _k |
| Immunogen: | Synthetic peptide SGGVVIA |
| Specificity: | β -amyloid 1-42 |
| Purification: | purified by affinity chromatography |
| Purity: | > 95% by SDS-PAGE |
| Solution: | PBS pH 7.4 without additives |
| Storage: | at -20 °C (repeated thawing and freezing should be avoided) |
| Applications: | ELISA, Western Blot, Immunohistochemistry |
| Dilution: | 2 ng/ml in ELISA, 1 μ g/ml in Western Blot |
| Reactivity: | β -amyloid(42), specific to C-terminal amino acids 41+42 |



Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

CEO
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig