

Anti-human phospho-231 TAU mab 8G5

Catalogue number	Package size
0102008701	100 µg
0102008703	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	8G5
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Synthetic peptide KKVAVVRT(Pi)PPKSPSS
Specificity:	human TAU phosphorylated at amino acid position 231 (T)
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot, Immunohistochemistry
Dilution:	1 µg/ml in Western Blot and IHC
Reactivity:	phosphorylated TAU at amino acid position 231 (T)

