

Anti-human TAU-E2/E3 mab 9E11

Catalogue number

847-0102006601
847-0102006603

Package size

100 µg
1 mg

For research use only.

Product:	monoclonal antibody
Clone:	9E11
Host:	mouse
Isotype:	IgGM
Immunogen:	Recombinant human TAU 441 peptides
Specificity:	Exon2/Exon 3 of human TAU441
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	Phosphate buffered saline ph 7.4, without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot
Dilution:	1 µg/ml in Western Blot
Reactivity:	Tau441, Tau410

