Revision: 003 Date: 28.10.2026

Immuno Beads Human p50-Tau (hTau-p50)

Catalogue number	Package size
847-0601000201	6 x IP
847-0601000202	24 x IP

For research use only.

Product:	Magnetic beads coupled antibody
Isotype:	Mouse IgG
Specificity:	Tau protein phosphorylated at amino acid T50
Solution:	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
Concentration:	10 mg magnetic beads/ml and 200 µg IgG/ml
Storage:	2-10°C *After delivery of the beads, make sure that the beads are completely in suspension and do not stuck to the tube walls or in the lid. To do this, vortex the beads and shortly spin them down by centrifuge so that all beads and storage liquid are at the bottom of the tube.*
Use for:	Immunoprecipitation of human p50-Tau protein
Application:	25 µl Immuno Beads per sample
Use with:	NEURO-IP Kit (847-0108000108)

Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

CEO
Dr. Ingolf Lachmann

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig