Revision: 002 Date: 28.10.2026

Immuno Beads Human patho-oligomeric α-Synuclein (hSCNA patho)

Catalogue number	Package size
847-0601000501	6 x IP
847-0601000502	24 x IP

For research use only.

Product:	Magnetic beads coupled antibody
Isotype:	Mouse IgG
Specificity:	Human α-Synuclein oligomers only, amino acids 47-52 (EGVVHGVA) in β-sheet
Solution:	PBS ph 7.4 containing BSA, Tween 20 and Na-Azide (0.02 %)
Concentration:	10 mg magnetic beads/ml and 200 µg IgG/ml
Storage:	2-10°C *After delivery of the beads, make sure that the beads are completely in suspension and do not stuck to the tube walls or in the lid. To do this, vortex the beads and shortly spin them down by centrifuge so that all beads and storage liquid are at the bottom of the tube.*
Use for:	Immunoprecipitation of human α-Synuclein oligomers
Application:	25 µl Immuno Beads per sample
Use with:	NEURO-IP Kit (847-0108000108)

Roboscreen GmbH
Hohmannstraße 7
04129 Leipzig • Germany

Phone: +49 341 989 734 0
Fax: +49 341 989 734 199
info@roboscreen.com
www.roboscreen.com

CEO
Dr. Ingolf Lachmann

Ust-IdNr. DE 813064314
County Court Leipzig
HRB 17503
Leipzig